




CONCENTRATORS CATALOG



CONCENTRATORS CATALOG

Concentrator CA-VH5

Concentrator CA-HG



Concentrator CA-V

Concentrator CA-PHT10

CONCENTRATOR CA-VH5



Model	CA-VH5
Flow rate (L/min)	0-5 L/min
Power consumption (W)	400W
Operating voltage	230V 50/60Hz
Concentration (%)	>90%
Outlet pressure (Mpa)	0.04-0.08 MPa
Sound level (dB)	<48 dB
Electrical category	Class II Typ BF
Net weight (kg)	17 kg
Dimension (mm)	W375xD252xH574 mm
Power failure alarm	When the machine switches off, a visual alarm is triggered
Temperature alarm	When the machine has more than 46 °C, the alarm is triggered
Pressure fault alarm	When the air pressure is low, the alarm is triggered
Spraying	Particle atomization <5μ reaches 90%; output pressure >0.1 Mpa
Low purity alarm	When the purity is <82% (+3), an alarm is triggered

CONCENTRATOR

CA-HG



Model	CA-HG 5	CA-HG 10
Power consumption (W)	390	480
Operating voltage	230V 50/60Hz	
Flow rate (L/min)	0-5 L/min	0-10 L/min
Concentration (%)	>90	
Outlet pressure (Mpa)	0.04-0.08	
Sound level (dB)	<45	<50
Electrical category	Class II Typ BF	
Product category	Class II b	
Net weight (kg)	27	30
Dimension (mm)	W400xD365xH650	
IP classification	IP 21	
Aerosol mask function (optional function)	<5 reaches 90% for spray type only	
Oxygen mask function (optional function)	When the oxygen purity is lower than 82%, an alarm is triggered (only for the oxygen purity alarm type)	

CONCENTRATOR CA-V



Model	CA-V8
Power consumption (W)	650
Operating voltage	230V 50/60Hz
Flow rate (L/min)	0-8 L/min
Concentration (%)	>90
Outlet pressure (Mpa)	0.04-0.08
Sound level (dB)	<50
Electrical category	Class II Typ BF
Product category	Class II a
Net weight (kg)	24
Dimension (mm)	W390xD337xH620
IP classification	IP 21
Aerosol mask function (optional function)	<5 reaches 90% for spray type only
Oxygen mask function (optional function)	When the oxygen purity is lower than 82%, an alarm is triggered (only for the oxygen purity alarm type)

CONCENTRATOR CA-PHT10



Name	portable oxygen concentrator
Model	CA-HPT10
Pulse transmission	transmission 1-5
Concentration	>90% (after 5 minutes of running)
Sound level (dB)	<58dB (A)
Decibel alarm	>60dB (A)
Nominal oxygen capacity	1-5 L/min
AC power supply	AC: 220V 50Hz
Battery	DC: 14.8V 6.4Ah 95Wh
Battery life	In gear 1, it is used continuously for 4 hours; In gear 5, continuously for 1.5 hours
Battery charging time	No more than 4 hours



RIMAVSKÁ SOBOTA
Malohontská 2594
+421 907 663 257
info@oxymedic.sk
oxymedic.sk